

081105T4DAP

DAIRY PLANT MANAGEMENT LEVEL 5

DA/OS/PM/CR/01/5/B

Procure Raw Milk

Nov/Dec 2024



**TVET CURRICULUM DEVELOPMENT, ASSESSMENT AND CERTIFICATION
COUNCIL (TVET CDACC)**

PRACTICAL ASSESSMENT

Time: 2 HOURS

INSTRUCTIONS TO CANDIDATE:

1. You are required to perform the following tasks
 - i. Prepare resazurin solution.
 - ii. Carry out Resazurin test on collected milk.
2. You have been provided with the following resources for the practical tasks:
 - Raw milk 100ml.
 - 1 Resazurin tablet
 - Distilled water
 - Measuring cylinder 100 ml
 - Hand Towel.
 - Thermometer

- Test tubes
- Pipette.
- Basin
- Milk jar
- Resazurin test interpretation chart
- Water
- Source of heat

easyvet.com